


micrux
TECHNOLOGIES




All-in-One platform for
thin-film electrodes



The innovative All-in-One cell (Ref. ED-AIO-CELL) provides an unique multipurpose interface with movable add-ons that can be easily interchanged for using the thin-film (micro)electrodes supplied by MicruX Technologies.

» AIO-cell general features

The AIO-cell enables the use of the thin-film (micro)electrodes in static (Drop / Batch-cell) or dynamic (Flow-cell) conditions, fulfilling the requirements of multiple electroanalytical applications.





The All-in-One cell is provided in two basic versions, the *base cell + one add-on* (Ref. ED-AIO-CELL-2x) or the *base cell + two add-ons* (Ref. ED-AIO-CELL-3x). Additional add-ons can be supplied with the base cell or by separated and they can be customized (design & materials) taking into account the specific requirements of the customers.

» AIO-cell specific Add-ons features

» Drop-cell (base cell)



The drop-cell provides a true user-friendly and robust interface with the instrumentation, enabling the use of microvolume (1 – 10 μL sample drops) with all standard (10 x 6 mm) thin-film (micro)electrodes. The drop-cell includes the basic assembly parts for the different add-ons.

» Batch-cell Add-ons



The batch-cell add-on enables the use of the standard (10 x 6 mm) thin-film (micro)electrodes in applications (batch analysis, standard additions, etc...) in which are required larger sample volume up to 400 μL .



» Flow-cell Add-ons



The flow-cell add-ons enable the use of standard (10 x 6 mm) thin-film (micro)electrodes as EC detection system in flowing liquids such as FIA, HPLC, CE, microfluidics, etc...



- » Wall-jet (inlet flow perpendicular to the working electrode surface) or Thin-layer (for microfluidic sensors) based flow-cell.
- » Standard fluidic ports ($\frac{1}{4}$ " - 28 UNF) with inlet channel of 0.5 mm I.D.
- » Low dead-volume (internal volume <500 nL or <60 nL). The cell volume is limited by an O-ring (2 or 3,5 mm I.D.) or microfluidic channels on chip.
- » High sensitivity electrochemical measurements.
- » Low sample requirements (microvolume < 20 μL).



Different **Flow-cell** and **Batch-cell add-ons** in **PMMA** (standard) and **PEEK** (on demand) are currently available for using in combination with the AIO platform and all standard (10 x 6 mm) thin-film (micro)electrodes.

» Methacrylate Add-ons

Standard add-ons are manufactured in methacrylate. PMMA is a suitable material for most of the analytical applications.



» Batch-cell Add-ons

- ❑ *Ref. BC-PMMA-2,0.* Compatible with thin-film (micro)electrodes with 2 mm diameter electrochemical cell. The cell is sealed with a 3 mm I.D. O-ring.
- ❑ *Ref. BC-PMMA-3,5.* Compatible with thin-film 10µElectrodes with 3,5 mm diameter electrochemical cell. The cell is sealed with a 3,5 mm I.D. O-ring.



» Flow-cell Add-ons (wall-jet & thin-layer)

- ❑ *Ref. FC-PMMA-2,0.* Compatible with thin-film (micro)electrodes with 2 mm diameter electrochemical cell. The cell volume is limited by a 2 mm I.D. O-ring.
- ❑ *Ref. FC-PMMA-3,5.* Compatible with thin-film 10µElectrodes with 3,5 mm diameter electrochemical cell. The cell volume is limited by a 3,5 mm I.D. O-ring.
- ❑ *Ref. TL-PMMA-1,0.* Compatible with microfluidic electrochemical sensors. The inlet/outlet of the microfluidic channel is sealed with two 1,15 mm I.D. O-rings.

» PEEK Add-ons

Add-ons are also available in PEEK (polyether ether ketone) on demand. PEEK offers advantages for applications where high mechanical and chemical resistance is required.



» Batch-cell Add-ons

- ❑ *Ref. BC-PEEK-2,0.* Compatible with thin-film (micro)electrodes with 2 mm diameter electrochemical cell. The cell is sealed with a 3 mm I.D. O-ring.
- ❑ *Ref. BC-PEEK-3,5.* Compatible with thin-film 10µElectrodes with 3,5 mm diameter electrochemical cell. The cell is sealed with a 3,5 mm I.D. O-ring.




» Flow-cell Add-ons (wall-jet & thin-layer)

- ❑ *Ref. FC-PEEK-2,0.* Compatible with thin-film (micro)electrodes with 2 mm diameter electrochemical cell. The cell volume is limited by a 2 mm I.D. O-ring.
- ❑ *Ref. FC-PEEK-3,5.* Compatible with thin-film 10µElectrodes with 3,5 mm diameter electrochemical cell. The cell volume is limited by a 3,5 mm I.D. O-ring.
- ❑ *Ref. TL-PEEK-1,0.* Compatible with microfluidic electrochemical sensors. The inlet/outlet of the microfluidic channel is sealed with two 1,15 mm I.D. O-rings.




» AIO-cell instrument-interfacing

» AIO-cell interfacing



Integrated easy-to-use connector (miniUSB-type connector) with pogo-pins

» Universal cable




The All-in-One cell is supplied with an universal cable compatible with any commercial electrochemical instrument.



miniUSB to AIO-cell

Female/male banana plug to instrument



Plug to instrumentation may be available in other format under previous request.



Severo Ochoa Building · Floor -1 – Room 4 & 6
Julián Clavería s/n · 33006 · Oviedo (Asturias) · SPAIN

Phone/FAX: +34 984151019

E-mail: info@micruxfluidic.com

Web: www.micruxfluidic.com

